

ACCOUNT & CLIENT INFORMATION	
Client Name	Date
Bill to # C50008642	Ship to
Costco	
Contact	Phone

To ensure proper fit, please indicate which RIC product this earmold will be used with:

KS8/KS7/Trax/Quintra/Legato/Adore/M-Core


FACEPLATE COLOR

Mocha (default) Beige Brown

SHELL COLORS	L	R
Tan (default)	<input type="checkbox"/>	<input type="checkbox"/>
Beige	<input type="checkbox"/>	<input type="checkbox"/>
Brown	<input type="checkbox"/>	<input type="checkbox"/>
Red	<input type="checkbox"/>	<input type="checkbox"/>
Blue	<input type="checkbox"/>	<input type="checkbox"/>
Clear	<input type="checkbox"/>	<input type="checkbox"/>

HP (75 DB) RECEIVER LENGTH	
L	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3
R	<input type="checkbox"/> 0 <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3

TIPS FOR MEASURING



- Choose the receiver length that corresponds to the top of the tragus.
- If in-between lengths, select the shorter length (larger number).

IF NEEDED, MAY WE CHANGE:

Venting	<input type="checkbox"/> Yes	<input type="checkbox"/> No	<input type="checkbox"/> Please Call
Canal Length	<input type="checkbox"/> Yes	<input type="checkbox"/> No	<input type="checkbox"/> Please Call

CANAL LENGTH

L	<input type="checkbox"/> Deep (default)	<input type="checkbox"/> Long	<input type="checkbox"/> Customer Select*
R	<input type="checkbox"/> Deep (default)	<input type="checkbox"/> Long	<input type="checkbox"/> Customer Select*

*Mark impression for proper length

SHELL OPTIONS/COATINGS	L	R
Clear Coat	<input type="checkbox"/>	<input type="checkbox"/>
Removal String	<input type="checkbox"/>	<input type="checkbox"/>
Flex Dip	<input type="checkbox"/>	<input type="checkbox"/>
Canal Grip	<input type="checkbox"/>	<input type="checkbox"/>
Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>

WAXGUARD	L	R
QuickGuard	<input type="checkbox"/>	<input type="checkbox"/>
Wax Spring ¹	<input type="checkbox"/>	<input type="checkbox"/>
Bell Bore ¹	<input type="checkbox"/>	<input type="checkbox"/>


¹QuickGuard cannot be combined with Bell Bore or Wax Spring

VENTING	L	R
No Vent	<input type="checkbox"/>	<input type="checkbox"/>
Micro Vent (1.0 mm) (default)	<input type="checkbox"/>	<input type="checkbox"/>
Pressure Vent (1.3 mm)	<input type="checkbox"/>	<input type="checkbox"/>
SAV Vent (1.6 mm)	<input type="checkbox"/>	<input type="checkbox"/>

SPECIAL INSTRUCTIONS

PLEASE MARK ANY AREA OF CONCERN

COMMENTS:



Canal Length Short Medium Long As Marked

Ear Texture Soft Medium Rigid NOTE: Please print clearly!