

APPLES TO APPLES

REXTON BICORE RUGGED BTE

RELY ON
LIFEPROOF
REXTON



SPEECH PRESERVATION TECHNOLOGY™		Rexton Rugged	Competitor
Soundpro 2.0	Automatically adjusts the amount of speech and noise separation to make listening easy in all situations	Yes	
Wireless Sync 4.0	Ultra-quick binaural processing preserves binaural cues even in dynamic environments to keep wearers hearing well in nearly all listening environments	Yes	
My Voice 2.0	Automatically reduces the loudness of your own voice when speaking	Yes	
AutoEcho Reducer	Improves comfort and speech understanding in reverberant places	Yes	
CONNECTIVITY			
HandsFree Mode	Allows wearer to talk on iPhones without holding the phone	Yes	
CallControl	Allows iPhone users to answer/end phone calls using the rocker switch on their hearing aid	Yes	
Rexton Assist*	Applies real-time improvements in hearing aid performance based on the wearer's needs	Yes	
Made for iPhone® & Android™ (ASHA)	Stream phone calls and audio from iPhone and compatible Android phones	Yes	
Streaming	Direct stream from your TV and other audio devices using Smart Transmitter 2.4 or Smart Mic	Yes	
SOUND QUALITY & COMFORT			
Multi-track Signal Processing	Follow conversations with ease, even in difficult listening situations	Yes	
Channels	Multi-track processing doubles the number of channels to 96 per hearing aid	Yes	
Motion Sensors	Automatic adjustment of directionality based on wearer movement	Yes	
Stereo iLock	Narrows the focus to hear a conversation partner directly in front of you when in noisy environments	Yes	
Auto iFocus 360	Microphone focus will automatically track to the direction of the loudest speech in a car environment	Yes	
Directional Hearing Control*	Allows you to manually customize the focus of the microphones	Yes	
Premium Noise Reduction	Reduces environmental noise such as fans, street noise, running water	Yes	
Intelligent Wind Noise Cancellation	Improves speech clarity when outdoors in windy environments and reduces wind noise	Yes	
Intelligent Feedback Preventer	Keeps hearing aids from whistling	Yes	
Sound Smoothing	Automatically reduces sudden noises, like dishes clanking	Yes	
Automatic Environment Manager	Automatically recognizes and adjusts to different environments including the car	Yes	
Bandwidth Compression	Allows better understanding of high pitch sounds for severe high frequency hearing losses	Yes	
Music Enhancer Programs	Includes optimized settings for musicians and listening to live or recorded music	3 programs	
XPhone	Allows you to hear phone conversation in both hearing aids while holding the phone up to one ear	Yes	
Telecoil/Loop	Improves hearing of speech delivered through a telecoil looped location and landline phones	Yes	
CONVENIENCE & ACCESSORIES			
Completely Automatic	Adapts to unique listening needs automatically as your environment changes	Yes	
Rechargeable	Lithium ion for up to 28 hours	Lithium ion	
Wearing Time	Including 5 hours streaming	24 hours	
Rexton App	Allows discreet control of hearing aid functions, remote tuning adjustments and Rexton Assist	Free	
Mask Mode*	Improves hearing voices muffled by face coverings and facemasks	Yes	
Programs	Allows the addition of customized programs if desired	6 programs	
Programmable Rocker Switch	On-board control of the hearing aid for volume and program change	Yes	
SecureTec Protection IP68	Keeps dust/dirt and moisture out of your hearing aids	Yes	
LifeProof™	Shower proof, 2x sweat and humidity resistance, drop/shock resistant	Yes	
Standard Charger	Easy access charger with open top	Included	
Smart Key	Discreet remote control to change volume and programs	Optional	
Smart Transmitter 2.4	Allows for direct streaming of Dolby Digital stereo sound from TVs	Optional	
Smart Mic	Remote microphone that allows streaming from Bluetooth enabled devices, offers hands-free mobile headset convenience and can be used as a remote control	Optional	
Price of 1 hearing aid, 1 charger, 3 year warranty, 2 year L&D			

*Available via the Rexton App

Rexton is a registered trademark of WSAUD A/S. SIV-16562/0423



Apple App Store is a trademark of Apple Inc.
iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.
Android and Google Play are trademarks of Google Inc.

HEARING.REXTON.COM

REXTON