

SUMMER 2023 SOLUTION ENHANCEMENTS

YOU ASKED, REXTON LISTENED



Based on your input we have made several improvements to your day-to-day workflow and the member experience.



New Length 5 Receiver

You requested a longer receiver to fit more members comfortably. The new length 5 MiniReceiver 3.0 is 4 mm longer than the L4 miniReceiver and is available in all receiver powers: S, M, P and HP (with custom shell).



Improved Aesthetics

The V2 miniReceiver 3.0 is 1 mm longer at the top curve of the receiver which results in the hearing aid sitting back a little further for a better fit. This will be a gradual change and these receivers can be identified by the V2 on the packaging.



New Concha Locks

We have returned to the sleeve style. Simply slip it over the receiver and attach the dome or sleeve. Available in S, M and P.



New Eco-friendly Packaging

As part of our shared commitment to the environment, we have changed our packaging to produce less waste and consume less cardboard. BiCore C R-Li and Li T bundles will be shipped in our current travel pouch covered by a stylish sleeve.



Phone Compatibility Check

Scan our QR code with your members' phones to quickly identify if they can stream calls and use the Rexton app.



Rexton Online Order Status

Rexton now offers the option for you to place orders online for hearing aids, accessories, supplies and repairs/remakes.

Your orders are confirmed immediately, ensuring order accuracy and improving delivery time. This allows us to get your orders to you faster.

Find it on Chrome under your Managed Bookmarks in the Vendor Order Status menu.

Rexton is a registered trademark of WSAUD A/S. SIV-16805 8/23.