

APPLES TO APPLES

REXTON REACH RECHARGEABLE RIC



SPEECH PRESERVATION TECHNOLOGY™		Rexton Reach	Competitor
Multi-Voice Focus	Allows you to focus and hear multiple speakers at the same time - even when there is competing background noise.	Yes	
Soundpro 3.0	Separates and adjusts the amount of speech and noise to make listening easy in all situations automatically	Yes	
Wireless Sync 4.0	Ultra-quick binaural processing ensures you are hearing well from all directions	Yes	
Voice Stabilizer	Helps you to keep hearing the voices you want when the environment around you gets louder	Yes	
My Voice 2.0	Automatically reduces the loudness of your own voice when speaking without affecting other sounds	Yes	
AutoEcho Reducer	Improves comfort and speech understanding in reverberant places	Yes	
CONNECTIVITY			
Made for iPhone® & Android™ (ASHA)	Stream phone calls and audio from iPhone and ASHA compatible Android phones	Yes	
Bluetooth® LE Audio Ready*	Will connect with mobile phones and other audio devices that support BLE Audio (firmware update required)	Yes	
Headset Mode*	Allows you to talk on compatible iPhones without holding the phone	Yes	
CallControl*	Allows iPhone users to answer/end phone calls using the rocker switch on their hearing aid	Yes	
Rexton Assist**	Allows you to change your hearing aid programming to improve sound quality right when you need it	In App	
Streaming	Direct stream from your TV and other audio devices using Smart Transmitter 2.4 or Smart Mic	Yes	
Remote Services**	Contact your Hearing Care Professional and receive programming assistance via Chat or a live audio/video connection	In App	
SOUND QUALITY & COMFORT			
Multi-track Signal Processing	Follow conversations with ease, even in difficult listening situations	Yes	
Channels	Multi-track processing doubles the number of channels to 96 per hearing aid	Yes	
Motion Sensors	Automatic adjustment of directionality based on wearer movement	Yes	
Stereo iLock	Narrows the focus to hear a conversation partner directly in front of you when in noisy environments	Yes	
Auto iFocus 360	Microphone focus will automatically track to the direction of the loudest speech in a car environment	Yes	
Directional Hearing Control**	Allows you to manually customize the focus of the microphones	In App	
Premium Noise Reduction	Reduces environmental noise such as fans, street noise, running water	Yes	
Wind Noise Canceller	Improves speech clarity when outdoors in windy environments and reduces wind noise	Yes	
Intelligent Feedback Preventer	Keeps hearing aids from whistling	Yes	
Sound Smoothing	Automatically reduces sudden noises, like dishes clanking	Yes	
Automatic Environment Manager	Automatically recognizes and adjusts to different environments including the car	Yes	
Bandwidth Compression	Allows better understanding of high pitch sounds for severe high frequency hearing losses	Yes	
Music Enhancer Programs	Includes optimized settings for musicians and listening to live or recorded music	3 programs	
Xphone	Allows you to hear phone conversation in both hearing aids when any phone is held up to one ear	Yes	
Telecoil/Loop	Improves hearing of speech delivered through a telecoil looped location and landline phones	Yes	
CONVENIENCE & ACCESSORIES			
Completely Automatic	Adapts to unique listening needs automatically as your environment changes	Yes	
Rechargeable	Lithium ion for up to 39 hours of runtime	Lithium ion	
Wearing Time	Including 5 hours streaming	34 hours	
Rexton App	Allows discreet control of hearing aid functions, remote tuning adjustments and Rexton Assist	Free	
Mask Mode**	Improves hearing voices muffled by face coverings and facemasks	In App	
Programs	Allows the addition of customized programs if desired	6 programs	
Programmable Rocker Switch	On-board control of the hearing aids for volume and program change	Yes	
ATEX Certified	Certified to be intrinsically safe to wear in explosive environments	Yes	
SecureTec Protection IP68	Keeps dust/dirt and moisture out of your hearing aids	Yes	
Travel Charger	Allows 3 full charges without needing to connect to a power outlet	Included	
Charging+ Station	Charges, dries and has UV sanitizing cycle	Optional	
Smart Transmitter 2.4	Allows for direct streaming of Dolby Digital stereo sound from TVs	Optional	
Smart Mic	Remote microphone allows streaming from Bluetooth devices, offers hands-free mobile headset convenience and can be used as a remote control	Optional	
Price of 2 hearing aids, travel charger, 3 year warranty, 2 year L&D			

*Requires iOS 15.2 and iPhone 11 or higher.

**Requires Rexton App



HEARING.REXTON.COM

Apple App Store is a trademark of Apple Inc.
 iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.
 Android and Google Play are trademarks of Google Inc.

REXTON

Rexton is a registered trademark of WSAUD A/S. SIV-19088-24