


APPLES TO APPLES

M-CORE LI-ION RECHARGEABLE BTEs



PLATFORM	M-Core BTEs			COMPETITOR
	B-Li M	B-Li P	B-Li HP	
Direct Streaming with Apple* iPhone®, iPad®, iPod®	Yes	Yes	Yes	
Bluetooth Ready for Android™ (with Smart Mic)	Yes	Yes	Yes	
Integrated Motion Sensor	Yes	Yes	Yes	
Processing Channels	48	48	48	
Completely Automatic	Yes	Yes	Yes	
Battery	Lithium-ion	Lithium-ion	Lithium-ion	
Wearing Time (No Streaming)	24	30	61	
Wearing Time (5 hours streaming)	21	27	57	
Programs	6	6	6	
Output/gain ThinTube	125/60	130/66	-	
Output/gain Earhook	133/60	135/77	140/82	
AUDIO EXCHANGE				
Auto Stereo iLock	Yes	Yes	Yes	
Auto iFocus 360/iFocus 360	Yes	Yes	Yes	
Directional Hearing Control (with Rexton App)	Yes	Yes	Yes	
XPhone	Yes	Yes	Yes	
SOUND QUALITY				
My Voice	Yes	Yes	Yes	
Voice Ranger	Yes	Yes	Yes	
Dynamic Extender	Yes	Yes	Yes	
Reverb Reducer Program	Yes	Yes	Yes	
Music Enhancer Programs	3	3	3	
TeleCoil/Loop Programs	-	Yes	Yes	
Premium Noise Reduction	Yes	Yes	Yes	
Intelligent Wind Noise Cancellation	Yes	Yes	Yes	
Intelligent Feedback Preventer	Yes	Yes	Yes	
Sound Smoothing	Yes	Yes	Yes	
CONVENIENCE				
Charging Station	Included	Included	Included	
Charging+ Station (charges, dries + UV sanitation)	-	Optional	Optional	
SecureTec Protection IP68	Yes	Yes	Yes	
Programmable Rocker Switch	Yes	Yes	Yes	
Volume Control and Program Change	Yes	Yes	Yes	
Smartphone App (Rexton App)	Free	Free	Free	
Remote Control (Smart Key)	Optional	Optional	Optional	
TV Transmitter (TV Sound)	Optional	Optional	Optional	
Remote Mic (Smart Mic)	Optional	Optional	Optional	
TeleCare Ready with Rexton App	Yes	Yes	Yes	
Price of 1 hearing aid, 3 year warranty, 2 year L&D				

Rexton is a registered trademark of WSAUD A/S. RX-22402



HEARING.REXTON.COM

*Apple App Store is a trademark of Apple Inc.
iOS is a trademark or registered trademark of Cisco in the U.S.
and other countries and is used under license.
Android and Google Play are trademarks of Google Inc.

REXTON