

## HEARING SOLUTIONS

REACH FEATURES	Reach Models		
	R PLUS RIC	RUGGED BTE	inoX-CIC
Multi-Voice Focus	✓	✓	
Multi-Track Processing	✓	✓	
Rechargeable	✓	✓	✓
Stream calls and audio from compatible Apple® and Android™ devices	✓	✓	
Hands-free/Tap Control	✓	iOS* Android**	
TeleCoil	✓	✓	
Bluetooth® Classic	✓		
MFi / ASHA / Bluetooth LE Audio		✓	
Mild to moderate hearing loss	✓	✓	✓
Mild to severe hearing loss	✓	✓	

\* Requires iPhone 11 and OS 15.2 or newer.  
\*\* With LE Audio compatible Android phones.

## ACCESSORIES\*



Our **Travel Chargers** are so compact you can take them anywhere, but powerful enough to provide 3-4 days of use on a single charge.



The **Charging+ Station** has an integrated drying function and UV cleaning light for extra care and protection.



**TV Sound** delivers high-quality audio directly from your TV to your Rexton hearing aids.



**Smart Key** — A small, convenient remote control.

\* Accessories not compatible with all models.

## ABOUT REXTON

Since 1955 Rexton has been a dedicated hearing aid brand delivering proven technology you can depend on. At work, at home, no matter the weather or your environment, you can always **RELY ON REXTON.**



Distributor:  
WS Audiology USA, Inc.  
101 Wood Avenue South  
7th Floor  
Iselin, NJ 08830

Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.  
Android, Google Play and the Google Play logo are trademarks of Google Inc.

Rexton is a registered trademark of WSAUD A/S. PN 1024-8771 // RX-23937

# REXTON



**RELY ON  
REXTON IN  
ANY SITUATION**

## REXTON REACH:

Make sure every voice is heard

When multiple people are talking, voices drop in and out and move around—making it challenging for hearing aid wearers to catch every word. **Multi-Voice Focus** changes that by scanning your environment 1,000 times per second with 4 focus beams, automatically adapting to shifts in speaker positions and volumes. Built on our **Multi-Track Processing**, which separates speech and background sounds into two independent tracks, Rexton Reach hearing aids ensure speech stands out clearly from the noise. **My Voice 2.0** keeps your own voice sounding natural and comfortable, while **Auto Echo Reducer** makes reverberant spaces like kitchens much easier to navigate. With Rexton, you're part of every conversation again.



## REXTON REACH R PLUS IS AVAILABLE IN 4 COLORS



Black

Sandy Brown

Rose Gold

Dark Champagne

## REASONS TO RELY ON REXTON

### UNIVERSAL CONNECTIVITY FOR EVERYONE

Connect to your smartphone or favorite devices — stream shows, music, and enjoy hands-free calling with ease. Rexton offers one of the most comprehensive connectivity portfolios available, so there's a style perfectly suited to your needs.

### POWER THAT KEEPS YOU GOING

Enjoy industry-leading battery runtime so you can rely on your hearing aids from morning to night. And when you're on the move, the included travel charger with built-in power bank keeps you powered wherever you go.

### AI-POWERED PERSONALIZATION

Rexton Assist, our cloud-based AI assistant built into the Rexton App, provides intelligent 24/7 support to help you make tailored, real-time sound adjustments on the go. Simply answer a few quick questions, and Rexton Assist automatically personalizes your hearing aids' performance when and where you need it—no appointment necessary.

## DEPENDABLE SUPPORT

Visit your smartphone's app store and download the Rexton App.



Conveniently connect with your Costco Hearing Center professional.



### LIVE ASSIST APPOINTMENTS

In addition to in-person visits, you can schedule "Live-Assist" video appointments from the comfort of your own home.



### REMOTE FINE-TUNING ADJUSTMENTS

Request adjustments anytime to optimize your hearing aid settings remotely.

Check your mobile device compatibility using the QR code below.

